

HALO[®]

Hemorrhage Arresting Lever[®] Operated

TWO HANDED

TOURNIQUET

ONE HANDED



Watch Video



1. Detach Red Buckle from base, position strap high on injured limb. re-hook Red Buckle to base.

1. Leave the Red Buckle attached to base, slide strap high on the injured limb.

Less initial tensioning is required for an arm than for a leg. Muscular thighs may require more initial Main Strap tensioning.

DO NOT cover the lever with the Main Strap.


Main Strap must be firmly in place before using Lever.

Finger Pull Strap may be used to secure the end of the Main Strap if movement over rough terrain is required.

If bleeding has not been controlled, release the Lever first by peeling the Locking Strap, tighten the Main Strap slightly and re-apply Lever. If bleeding continues, apply a second HALO next to the first.




2. Free Main Strap so that it slides freely through Red Buckle.



3. Pull Main Strap away from base till snug. Use finger pull as necessary.



4. Wrap Main Strap back over itself around limb and re-attach by pressing down firmly.



5. Pull Locking Strap outward from lever. Using the palm, Push lever over, flat against the base.

6. Secure Lever with Locking Strap by pressing strap firmly onto the Main Strap.

7. Check for distal pulse and/or control of Hemorrhage.

8. Record Time on Time Strap.